

ABD Trauma Bandage

CE REF TM-212

NSN: 6510-01-541-8121

Vol. 40 per case, Total Weight: 7kg/15.4lbs

TyTek Medical's Abdominal Trauma Bandage is a superior pressure dressing providing haemorrhage control for abdominal and large extremity wounds in pre-hospital emergency care. The integrated large 12" x 12" non-stick pad makes this pressure dressing an ideal solution for managing massive wounds. It delivers a one-step fast application of the pressure dressing, eliminating the need for separate gauze pads and elastic wraps.

Unlike traditional elastic wraps which roll away or become uncontrollable in the field, our trauma bandages have integrated hook and loop stays throughout the bandage, ensuring controllability in all applications.


Details:

- Sterile, oversized 12" x 12" non-stick pad suitable for applying immediate pressure
- Easy-open white tear notches for quick access
- Intermittent hook and loop prevents pressure dressing from unrolling during application
- Simple, quick, controlled and effective pressure dressing system
- Vacuum sealed to minimize cubic space
- Most cost effective large wound trauma bandage on the market

Dimensions:

- Bandage: L 70" x W 6"
- Packaged: L 5 1/4" x W 6" x D 1 1/2"





Product Range Overview

Product Name	Part No.	NSN	CE Mark	Volume	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	100 per case	8kg/17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	50 per case	6kg/13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	40 per case	7kg/15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	200 per case	9kg/19.8lbs
TPAK	TM-303	6515-01-541-0635	Yes	1000 per case	15.5kg/34lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	800 per case	16kg/35.2lbs
Xtricate	TM-511	6530-01-477-8525	Yes	1 per case	8kg/17.6lbs
TVAC	TM-909	6515-01-540-7206	Yes	20 per case	4.8kg/10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	20 per case	1.6kg/3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	30 per case	1.4kg/3lbs