



## Piranha Trauma Shears Large (7.25")

CE  TM-707

NSN: 6515-01-538-9276


Vol. 20 per case, Total Weight: 1.6kg/3.5lbs


TyTek Medical's 7 1/4" Piranha Trauma Shears are the very best tactical shears available anywhere. These tactical shears were engineered to meet and exceed the military's stringent standards. Our trauma shears will cut through the toughest materials - web gear, denim, leather, boots, heavy winter clothing and all the metal components that come with them.

More importantly, TyTek Medical's tactical trauma shears are versatile - they easily cut the softest materials like gauze, which is where traditional trauma shears fail. Totally unique to TyTek Medical, is the ergonomic comfort grip handle with an inlayed rubber grip. It provides maximum comfort and control in high stress pre-hospital emergency care scenarios.

### Details:

- Premium-grade, durable stainless steel tactical shears blades for precision sharpness
- Reinforced, serrated edge for maximum cutting ability through the toughest materials
- Clothing-lift lip and rounded tip for patient safety
- Ergonomic inlayed rubber grip and contoured handles for maximum comfort and control





## Product Range Overview

Product Name	Part No.	NSN	CE Mark	Volume	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	100 per case	8kg/17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	50 per case	6kg/13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	40 per case	7kg/15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	200 per case	9kg/19.8lbs
TPAK	TM-303	6515-01-541-0635	Yes	1000 per case	15.5kg/34lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	800 per case	16kg/35.2lbs
Xtricate	TM-511	6530-01-477-8525	Yes	1 per case	8kg/17.6lbs
TVAC	TM-909	6515-01-540-7206	Yes	20 per case	4.8kg/10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	20 per case	1.6kg/3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	30 per case	1.4kg/3lbs