

## TPAK

CE  TM-303

NSN: 6515-01-541-0635


Vol. 1000 per case, Total Weight: 15.5kg/34lbs

TyTek Medical's 14 gauge TPAK chest decompression needle is in the recommended 3.25" length for treating a tension pneumothorax. You can't count on a standard IV catheter to needle a chest deeply enough in needle thoracostomy. Our 14 gauge x 3.25" TPAK chest decompression needle is a compact, reliable solution for treating a tension pneumothorax. Shock resistant case ensures the needle/catheter is secure upon deployment, and the capless flash chamber gives instant confirmation of needle placement in pre-hospital needle thoracostomy. TPAK comes in a rugged case, with a textured, easy-open hex cap.

### Details:

- 14 gauge x 3.25" chest decompression needle
- Strong, reliable tension pneumothorax needle allows for secure placement of catheter for continuous relief during needle thoracostomy
- Rugged tube with textured, easy-open hex cap
- Capless flash chamber for instant confirmation of needle placement
- Convenient, compact size for maximum portability
- Shock resistant, precision dimensioned case and needle hub ensure that chest decompression needle and catheter secure upon deployment for safety





## Product Range Overview

| Product Name                           | Part No. | NSN              | CE Mark | Volume        | Weight        |
|--|----------|------------------|---------|---------------|---------------|
| 4"Trauma Bandage                       | TM-204   | 6510-01-558-4108 | Yes     | 100 per case  | 8kg/17.6lbs   |
| 6"Trauma Bandage                       | TM-206   | 6510-01-558-4114 | Yes     | 50 per case   | 6kg/13.2lbs   |
| ABD Trauma Bandage                     | TM-212   | 6510-01-541-8121 | Yes     | 40 per case   | 7kg/15.4lbs   |
| EZ Gauze                               | TM-215   | 6510-01-605-2234 | Yes     | 200 per case  | 9kg/19.8lbs   |
| TPAK                                   | TM-303   | 6515-01-541-0635 | Yes     | 1000 per case | 15.5kg/34lbs  |
| PneumoDart                             | TM-317   | 6515-01-655-9514 | Yes     | 800 per case  | 16kg/35.2lbs  |
| Xtricate                               | TM-511   | 6530-01-477-8525 | Yes     | 1 per case    | 8kg/17.6lbs   |
| TVAC                                   | TM-909   | 6515-01-540-7206 | Yes     | 20 per case   | 4.8kg/10.5lbs |
| Piranha Trauma Shears<br>Large (7.25") | TM-707   | 6515-01-538-9276 | Yes     | 20 per case   | 1.6kg/3.5lbs  |
| Piranha Trauma Shears<br>Small (5.50") | TM-705   | 6515-01-538-9276 | Yes     | 30 per case   | 1.4kg/3lbs    |