

EZ Gauze

CE  TM-215

NSN: 6510-01-605-2234

Vol. 200 per case, Total Weight: 9kg/19.8lbs

Our compressed gauze is superior to typical rolled gauze. The unique folding and dispensing system makes our compressed gauze the ideal solution for fast and easy wound packing. The system allows our gauze to be easily used and cut with virtually no wasted gauze or time. In pre-hospital emergency care scenarios, rolled gauze can become damaged and lose its sterility.


TyTek Medical's compressed gauze maintains sterility in rugged, compact packaging that facilitates wound packing and wound dressing. In wound packing situations, gauze removal from the wound tract is made easier because EZ-Gauze is a continuous-feed gauze. It is removed in a single continuous piece, as opposed to removing an unknown number of 4x4 gauze pads from the wound.

Additionally, instead of separate clips or tabs, which can be easily lost, our Abdominal Trauma Bandage has integrated clips at the end of the bandage to easily secure the pressure dressing to the patient.

Details:

- 4 yards, 6-ply, 4.5" sterile compressed gauze
- Vacuum sealed for maximum portability
- Easy-open white tear notches for quick and easy access for wound packing and wound dressing
- Unique z-rolled dispensing system resulting in controlled application with no wasted gauze





Product Range Overview

Product Name	Part No.	NSN	CE Mark	Volume	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	100 per case	8kg/17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	50 per case	6kg/13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	40 per case	7kg/15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	200 per case	9kg/19.8lbs
TPAK	TM-303	6515-01-541-0635	Yes	1000 per case	15.5kg/34lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	800 per case	16kg/35.2lbs
Xtricate	TM-511	6530-01-477-8525	Yes	1 per case	8kg/17.6lbs
TVAC	TM-909	6515-01-540-7206	Yes	20 per case	4.8kg/10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	20 per case	1.6kg/3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	30 per case	1.4kg/3lbs