

## TPAK10

CE REF TM-310

NSN: 6515-01-673-1701


Vol. 1000 per case, Total Weight: 15.5kg/34lbs

TyTek Medical's 10 gauge TPAK chest decompression needle is in the recommended 3.25" length for treating a tension pneumothorax. You can't count on a standard IV catheter to needle a chest deeply enough in needle thoracostomy. Our 10 gauge x 3.25" TPAK chest decompression needle is a compact, reliable solution for treating a tension pneumothorax. Shock resistant case ensures the needle/catheter is secure upon deployment, and the capless flash chamber gives instant confirmation of needle placement in pre-hospital needle thoracostomy. TPAK comes in a rugged case, with a textured, easy-open hex cap.

### Details:

- 10 gauge x 3.25" chest decompression needle
- Strong, reliable tension pneumothorax needle allows for secure placement of catheter for continuous relief during needle thoracostomy
- Rugged tube with textured, easy-open hex cap
- Capless flash chamber for instant confirmation of needle placement
- Convenient, compact size for maximum portability
- Shock resistant, precision dimensioned case and needle hub ensure that chest decompression needle and catheter secure upon deployment for safety





## Product Range Overview

Product Name	Part No.	NSN	CE Mark	Volume	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	100 per case	8kg/17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	50 per case	6kg/13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	40 per case	7kg/15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	200 per case	9kg/19.8lbs
TPAK	TM-310	6515-01-541-0635	Yes	1000 per case	15.5kg/34lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	800 per case	16kg/35.2lbs
Xtricate	TM-511	6530-01-477-8525	Yes	1 per case	8kg/17.6lbs
TVAC	TM-909	6515-01-540-7206	Yes	20 per case	4.8kg/10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	20 per case	1.6kg/3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	30 per case	1.4kg/3lbs