

4" Trauma Bandage

CE  TM-204

NSN: 6510-01-558-4108

Vol. 100 per case, Total Weight: 8kg/17.6lbs

MOQ: Case Price: \$499.00

Our 4" Trauma Bandages are simple and effective pressure dressings providing hemorrhage control in pre-hospital emergency care situations. Our pressure dressings incorporate a non-stick pad affixed to the bandage, delivering a one-step fast application of the pressure dressing, eliminating the need for separate gauze pads and elastic wraps.

Unlike traditional elastic wraps which roll away or become uncontrollable in the field, our trauma bandages have integrated hook and loop that features throughout the bandage, ensuring controllability in all applications. Additionally, instead of separate clips or tabs, which can be easily lost, our pressure dressings have integrated clips at the end of the bandage to easily secure the trauma bandage to the patient.

Details:


- Sterile, non-stick pad suitable for applying immediate pressure to wounds
- Easy-open white tear notches for easy ID and quick access
- Intermittent hook and loop prevents trauma bandage from unrolling during application
- Simple, quick, controlled and effective pressure dressing system
- Vacuum sealed to minimize cubic space
- Most cost effective trauma bandage on the market

Dimensions:

- Bandage: L 70" x W 4"
- Packaged: L 4 1/2" x W 2 1/2" x D 1 3/4"



Rev 08262020



Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020