

ABD Trauma Bandage

CE  TM-212

NSN: 6510-01-541-8121

Vol. 40 per case, Total Weight: 7kg/15.4lbs

MOQ: Case Price: \$311.60

TyTek Medical's Abdominal Trauma Bandage is a pressure dressing providing hemorrhage control for abdominal and large extremity wounds in pre-hospital emergency care. The integrated large 12" x 12" non-stick pad makes this pressure dressing an ideal solution for managing massive wounds. It delivers a one-step fast application of the pressure dressing, eliminating the need for separate gauze pads and elastic wraps.

Unlike traditional elastic wraps which roll away or become uncontrollable in the field, our trauma bandages have integrated hook and loop that features throughout throughout the bandage, ensuring controllability in all applications.

Additionally, instead of separate clips or tabs, which can be easily lost, our Abdominal Trauma Bandage has integrated clips at the end of the bandage to easily secure the pressure dressing to the patient.

Details:


- Sterile, oversized 12" x 12" non-stick pad suitable for applying immediate pressure
- Easy-open white tear notches for quick access
- Intermittent hook and loop prevents pressure dressing from unrolling during application
- Simple, quick, controlled and effective pressure dressing system
- Vacuum sealed to minimize cubic space

Dimensions:

- Bandage: L 70" x W 6"
- Packaged: L 5 1/4" x W 6" x D 1 1/2"



Rev 08262020



Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020