

## Piranha Trauma Shears Small (5.50")

CE  TM-705

NSN: 6515-01-538-9276

Vol. 30 per box, Total Weight: 1.4kg/3lbs

MOQ: Box Price: \$260.70

TyTek Medical's 5 1/2" Piranha Trauma Shears are premium grade tactical shears. These tactical shears were engineered to meet and exceed the military's stringent standards. Our trauma shears will cut through the toughest materials - web gear, denim, leather, boots, heavy winter clothing and all the metal components that come with them.


More importantly, TyTek Medical's tactical trauma shears are versatile - they easily cut the softest materials like gauze, which is where traditional trauma shears fail. Totally unique to TyTek Medical, is the ergonomic comfort grip handle with an inlayed rubber grip. It provides maximum comfort and control in high stress pre-hospital emergency care scenarios.

### Details:

- Premium-grade, durable stainless steel tactical shears blades for precision sharpness
- Reinforced, serrated edge for maximum cutting ability through the toughest materials
- Clothing-lift lip and rounded tip for patient safety
- Ergonomic inlayed rubber grip and contoured handles for maximum comfort and control



Rev 08262020



## Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020