

TPAK 10g

CE  TM-310

NSN: 6515-01-673-1701

Vol. 1000 per case, Total Weight: 15.5kg/34lbs

MOQ: Box (50 Per Box) Price: \$499.50


TyTek Medical's 10 gauge TPAK chest decompression needle is the Committee on Tactical Combat Casualty Care's (CoTCCC) recommended 3.25" length and gauge for treating a tension pneumothorax.

Details:

- 10 gauge x 3.25" chest decompression needle
- Strong, reliable tension pneumothorax needle allows for secure placement of catheter for continuous relief during needle thoracostomy
- Rugged tube with textured, easy-open hex cap
- Capless flash chamber for instant confirmation of needle placement
- Convenient, compact size for maximum portability
- Shock resistant, precision dimensioned case and needle hub ensure that chest decompression needle and catheter secure upon deployment for safety



Rev 08262020



Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020