

## TPAK 14g

CE  TM-303

NSN: 6515-01-541-0635

Vol. 1000 per case, Total Weight: 15.5kg/34lbs

MOQ: Box (50 Per Box) Price: \$499.50


TyTek Medical's 14 gauge TPAK chest decompression needle is the Committee on Tactical Combat Casualty Care's (CoTCCC) recommended 3.25" length and gauge for treating a tension pneumothorax.

### Details:

- 14 gauge x 3.25" chest decompression needle
- Strong, reliable tension pneumothorax needle allows for secure placement of catheter for continuous relief during needle thoracostomy
- Rugged tube with textured, easy-open hex cap
- Capless flash chamber for instant confirmation of needle placement
- Convenient, compact size for maximum portability
- Shock resistant, precision dimensioned case and needle hub ensure that chest decompression needle and catheter secure upon deployment for safety



Rev 08262020



## Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020