

## TVAC

CE  TM-909

NSN: 6515-01-540-7206

Vol. 20 per case, Total Weight: 4.8kg/10.5lbs

MOQ: Case Price: \$399.80


TyTek Medical's disposable portable suction unit with its resilient silicone bulb performs reliably in even the harshest climates without loss of suction. The TVAC Manual Suction Unit comes in a field friendly compact package, just pull out the built-in tube and start sucking. No bulky cartridges or brittle parts to break. After use, simply place the entire TVAC Manual Suction Unit back inside its packaging, seal, expose the biohazard label on the package and dispose appropriately. This entirely removes the possibility of cross contamination.

### Details:

- Strong, silicone bulb creates vacuum pressure of approx. 110 mmHg
- Collection bag volume >1,000ml
- Resilient, reliable silicone bulb prevents loss of suction in even the harshest climates
- Contoured suction tube tip prevents tissue damage during application
- Reliable, simple and inexpensive manual suction unit for clearing airway obstructions
- Simple disposal eliminates cross contamination concerns



Rev 08262020



## Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020