



XTricate

CE  TM-511

NSN: 6530-01-477-8525

Vol. 1 each, Total Weight: 8kg/17.6lbs

MOQ: 1 Price: \$699.00

TyTek Medical's XTricate is our tech rescue extrication and immobilization sleeve. XTricate is ideal in confined space rescue and tech rescue scenarios. It is versatile in tech rescue situations, allowing both vertical and horizontal lifting and support. It will immobilize patients with neck or spinal injuries, as well as extricate patients in confined space rescue scenarios.


Extremely durable, XTricate has been independently lab tested to carry over 1,000lbs. One vertical lifting point and 4 horizontal lifting points for helicopter hoist capability. Five chest and leg straps and one fail safe strap are secured with Mil Spec buckles.

Details:

- Simple, secure tech rescue immobilization/extrication sleeve with 4 horizontal and 1 vertical Mil Spec lifting-strap attachments
- (1) Vertical lifting point D-Ring specs: 100% proof load tested @ 5,000lbs, Min tensile strength 5,000lbs, meets or exceeds PIA-H-7195, ANSI Z-359.2, and CSA Z-259.12-01
- (4) Horizontal lifting point O-Ring specs: 100% proof load tested @2,500lbs, meets or exceeds PIA-H-7195
- Payload: 1,000lbs
- Independent lab tested with vertical lift failure point of over 10,000lbs and tested to a horizontal failure point of over >7,400
- Spine board compartment
- Removable head blocks with medic easy-access holes
- Adjustable head and chin straps
- Zippered storage compartments
- 12 Padded carry handles
- Replacement parts available: Head blocks, chin strap, and carrying case
- Easy to clean and waterproof



Rev 08262020



Product Range Overview

Product Name	Part No.	NSN	CE Mark	MOQ	Price	Weight
4"Trauma Bandage	TM-204	6510-01-558-4108	Yes	Case (100 Per Case)	\$499.00	8kg / 17.6lbs
6"Trauma Bandage	TM-206	6510-01-558-4114	Yes	Case (50 Per Case)	\$284.50	6kg / 13.2lbs
ABD Trauma Bandage	TM-212	6510-01-541-8121	Yes	Case (40 Per Case)	\$311.60	7kg / 15.4lbs
EZ Gauze	TM-215	6510-01-605-2234	Yes	Case (200 Per Case)	\$918.00	9kg / 19.8lbs
TPAK 14ga	TM-303	6515-01-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.7lbs
TPAK 10ga	TM-310	6515-10-541-0635	Yes	Box (50 Per Box)	\$499.50	0.8kg / 1.8lbs
PneumoDart	TM-317	6515-01-655-9514	Yes	Box (40 Per Box)	\$1039.60	0.8kg / 1.7lbs
XTricate	TM-511	6530-01-477-8525	Yes	Each	\$699.00	8kg / 17.6 lbs
TVAC	TM-909	6515-01-540-7206	Yes	Case (20 Per Case)	\$399.80	4.8kg / 10.5lbs
Piranha Trauma Shears Large (7.25")	TM-707	6515-01-538-9276	Yes	Box (20 Per Box)	\$184.80	1.6kg / 3.5lbs
Piranha Trauma Shears Small (5.50")	TM-705	6515-01-538-9276	Yes	Box (30 Per Box)	\$260.70	1.4kg / 3lbs

Rev 08262020